







KATANA Zirconia (ML) Handy Chart

(ノリタケカタナジルコニア)

Thickness : Before Sintering (After Sintering : Approx.)


	14mm (11.2mm)	18mm (14.4mm)	22mm (17.6mm)
ML : A Light			
ML : A Dark			
ML : B Light			

Printed color can be slightly different from actual color.

KATANA Zirconia (ML : Multi Layered)


(ノリタケカタナジルコニア)

Specialized for FCZ (Full Contour Zirconia)



Sintering program

Max.Temp.1,500°C
(Hold time : 2hrs)



Kuraray Noritake Dental Inc.

300 Higashiyama, Miyoshi-cho
Miyoshi, Aichi 470-0293, Japan
<http://www.kuraraynoritake.com>

製造販売元

クラレノリタケデンタル株式会社

新潟県胎内市倉敷町2-28

ノリタケカタナジルコニア

管理医療機器 歯科切削加工用セラミックス
医療機器認証番号 223AFBZX00185000

印刷のため実際の色調と異なる場合があります。



CE 0120